


# FastSystem<sup>®</sup>

Upper Abdominal Retractor System

BF200

Table-Mounted Retraction for  
Unparalleled Upper Abdominal  
Exposure with Superior  
Subcostal Margin Elevation



## Excellent Retraction on Morbidly Obese Patients

- Two-post setup provides a rigid, low-profile frame
- Frame adjusts to fit operating room table extenders
- Optional obesity blade kit for increased visibility


## Fast Setup—Fast Mid-Surgery Changes

- Sets up in less than a minute
- Easy-to-use one piece components

# Product Specifications

## BF200 FastSystem® Upper Abdominal Retractor System

A	2	Omni-Clamp® Sterile Field Post, Double Joint	4944
B	1	Standard Cross Bar	4943
C	1	Hinged Cross Bar	4940
D	7	SnapClamp® Blade Holder	4040
E	2	Angled Side Arm, 8" x 18" (20.3cm x 45.7cm)	3759
F	1	Doherty-Harrington Retractor Blade, 5 1/2" x 3 1/2" (14cm x 8.9cm)	3091
G	2	Sternal Retractor, 3" x 3" (7.6cm x 7.6cm)	3094
H	2	Mayo Swivel Retractor, 3 1/2" x 2 3/4" (8.9cm x 7cm)	3536
I	1	Blade Holder, 9" (22.9cm)	3623
J	1	Malleable Swivel Retractor, 2" x 6 1/2" (5cm x 16.5cm)	3872
K	1	Malleable Swivel Retractor, 3" x 6 1/2" (7.6cm x 16.5cm)	3876



4849 White Bear Parkway, St. Paul, MN 55110 USA  
 800-367-8657 • 651-287-4300 • fax: 651-287-4400  
[www.omni-tract.com](http://www.omni-tract.com)

© 2006 Omni-Tract Surgical  
 All Rights Reserved.

U.S. Patent Nos. 4,355,631; 4,617,916; 4,718,151; 4,930,932; 4,949,707; 5,020,195;  
 5,242,240; 5,380,338; 5,400,772; 5,688,223; 5,727,899; 5,792,046; 5,882,298;  
 5,899,627; 6,264,396; 6,315,718; D343,235; Other Patents Pending.  
 Omni-Tract, FastSystem, Omni-Clamp, WISHBONE and SnapClamp  
 are registered trademarks of Minnesota Scientific, Inc.



1-06/40-6162/2500