


Obese Patients

### Obesity Blade Kit — OK200



*The Obesity Blade Kit has deeper retractors for use on large or obese patients*

OK200 Obesity Blade Kit			
A	2	Splanchnic Retractor, 2" x 8" (5 cm x 20.3 cm)	3529
B	2	Renal Vein Retractor, 1" x 8" (2.5 cm x 20.3 cm)	4929
C	1	Deaver Retractor, 2-3/4" x 8" (7 cm x 20.3 cm)	3531
D	2	Mayo Swivel Retractor, 3-1/2" x 2-3/4" (8.9 cm x 7 cm)	3536
E	1	Fence Swivel Retractor, 4" x 8" (10.2 cm x 20.3 cm)	3856
F	1	Malleable Swivel Retractor, 2" x 8" (5 cm x 20.3 cm)	3874
G	1	Malleable Swivel Retractor, 3" x 8" (7.6 cm x 20.3 cm)	3878
H	1	Malleable Swivel Retractor, 4" x 8" (10.2 cm x 20.3 cm)	3882

### Extension Kit — EK300

*The Extension Kit completely surrounds the operative site for more precise blade placement and extra frame rigidity*



EK300 Extension Kit			
A	1	Hinged Extension Arm	3588
B	2	SnapClamp	4040



4849 White Bear Parkway  
St. Paul, MN 55110 USA  
1-800-367-8657  
651-287-4300  
fax: 651-287-4400  
<http://www.omni-tract.com>

