


FastSystem®

Pelvic Retractor Systems



UF100

Table-Mounted Retraction for Hands-Free Exposure during Vaginal, Perineal and Minimally Invasive Abdominal Approaches

Easy Access to Conventionally-Limited Operative Field

- Small, independently-adjustable, low-profile frame arms follow contours of the operative field
- Multi-plane, multi-positional hands-free retraction
- Liberates surgical staff from awkward positioning when holding retractors


Fast Setup—Fast Mid-Surgery Changes

- Sets up in less than a minute
- Extension arm allows convenient system placement anywhere on table
- Open WISHBONE® Frame allows unobstructed drainage and use of weighted speculum/TUR drapes
- Easy-to-use one piece components

Product Specifications

UF100 FastSystem® Pelvic Retractor System

A	1	Omni-Clamp® Sterile Field Post	4080
B	1	FastSystem® Small WISHBONE® Frame	4985
C	5	SnapClamp® Blade Holders	4040
D	1	Extension Arm 20" (51cm)	4988
E	4	Malleable Swivel Retractor, 3/4" x 2" (1.9cm x 5cm)	4964
F	4	Malleable Swivel Retractor, 3/4" x 3" (1.9cm x 7.6cm)	4965
G	4	Malleable Swivel Retractor, 3/4" x 4" (1.9cm x 10.2cm)	4966
H	1	Posterior Swivel Retractor, 1 3/4" x 3" (4.5cm x 7.6cm)	4975
I	1	Posterior Swivel Retractor, 1 3/4" x 4" (4.5cm x 10.2cm)	4976
J	1	Anterior Swivel Retractor, 1 1/2" x 3" (3.8cm x 7.6cm)	4977



4849 White Bear Parkway, St. Paul, MN 55110 USA
 800-367-8657 • 651-287-4300 • fax: 651-287-4400
www.omni-tract.com

© 2005 Omni-Tract Surgical
 All Rights Reserved.

U.S. Patent Nos. 4,355,631; 4,617,916; 4,718,151; 4,930,932; 4,949,707; 5,020,195;
 5,242,240; 5,380,338; 5,400,772; 5,688,223; 5,727,899; 5,792,046; 5,882,298;
 5,899,627; 6,264,396; 6,315,718; D343,235; Other Patents Pending.
 Omni-Tract, FastSystem, Omni-Clamp, WISHBONE and SnapClamp
 are registered trademarks of Minnesota Scientific, Inc.



8-05/40-6140/2500