



Omni-Flex™

General Surgery Retractor System



Your Patient
Your Approach

Product Specifications



Omni-Flex™ General Retractor Systems

G0450 Expanded General Retractor System

G0400 General Retractor System

G0450	G0400	Description	Part Number
A	1	1 Sterile Field Post	10244
B	1	1 Support Arm	11050
C	2	2 Standard WISHBONE® Frame Arm	10269
D	--	-- Small WISHBONE® Frame Arm (Optional)	10270
E	8	6 SideSnap™ Clamp	11030
F	1	-- Hinged Extension Arm, 19" (48.3cm)	11121
G	1	1 Deaver Swivel Retractor, 2 3/4" x 4 1/2" (7cm x 11.4cm)	3784
H	1	1 Harrington Sweetheart Swivel Retractor, 2" x 5" (5cm x 12.7cm)	3794
I	1	1 Fence Swivel Retractor, 4" x 6" (10.2cm x 15.2cm)	3772
J	1	1 Kelly Modified Swivel Retractor, 3" x 2 1/2" (7.6cm x 6.4cm)	3862
K	1	1 Kelly Modified Swivel Retractor, 3" x 4" (7.6cm x 10.2cm)	3864
L	1	1 Malleable Slotted Swivel Retractor, 2" x 6 1/2" (5cm x 16.5cm)	3788
M	1	1 Malleable Slotted Swivel Retractor, 3" x 6 1/2" (7.6cm x 16.5cm)	3790
N	2	2 Mayo Swivel Retractor, 2 3/4" x 2" (7cm x 5cm)	3335
O	1	1 Splanchnic Swivel Retractor, 2" x 6" (5cm x 15.2cm)	3774



4849 White Bear Parkway, St. Paul, MN 55110 USA
 800-367-8657 • 651-287-4300 • fax: 651-287-4400
 www.omni-tract.com

© 2006 Omni-Tract Surgical
 All Rights Reserved.

U.S. Patent Nos. 4,355,631; 4,617,916; 4,718,151; 4,930,932; 4,949,707;
 5,020,195; 5,242,240; 5,380,338; 5,400,772; 5,688,223; 5,727,899;
 5,792,046; 5,882,298; 5,899,627; 6,264,396; 6,315,718; 6,663,563,
 D343,235; Other Patents Pending. Omni-Flex and SideSnap are trademarks;
 Omni-Tract and WISHBONE are registered trademarks of Minnesota
 Scientific, Inc.



5-06/40-6180/2500