



Omni-Flex™

Urology Retractor System



Your Patient
Your Approach

Product Specifications



Omni-Flex™ Urology Retractor Systems

U0450 Expanded Urology Retractor System
U0400 Urology Retractor System
U0100 Pelvic Retractor System

U0450	U0400	U0100	Description	Part Number
A	1	1	Sterile Field Post	10244
B	1	1	Support Arm	11050
C	2	2	Standard WISHBONE® Frame Arm	10269
D	--	--	Small WISHBONE® Frame Arm	10270
E	8	6	SideSnap™ Clamps	11030
F	1	--	Hinged Extension Arm, 19" (48.3cm)	11121
G	--	--	Extension Arm, 20" (50.8cm)	11070
H	--	--	Anterior Swivel Retractor, 1½" x 3" (3.8cm x 7.6cm)	4977
I	1	1	Deaver Swivel Retractor, 2¾" x 4½" (7cm x 11.4cm)	3784
J	1	1	Deaver Swivel Retractor, 1½" x 6" (3.8cm x 15.2cm)	3854
K	1	1	Fence Swivel Retractor, 4" x 6" (10.2cm x 15.2cm)	3772
L	1	1	Harrington Sweetheart Swivel Retractor, 2" x 5" (5cm x 12.7cm)	3794
M	1	1	Harrington Sweetheart Swivel Retractor, 2½" x 7" (6.4cm x 17.8cm)	3796
N	--	--	Malleable Swivel Retractor, ¾" x 2" (1.9cm x 5cm)	4964
O	--	--	Malleable Swivel Retractor, ¾" x 3" (1.9cm x 7.6cm)	4965
P	--	--	Malleable Swivel Retractor, ¾" x 4" (1.9cm x 10.2cm)	4966
Q	1	1	Malleable Slotted Swivel Retractor, 2" x 6½" (5cm x 16.5cm)	3788
R	1	1	Malleable Slotted Swivel Retractor, 3" x 6½" (7.6cm x 16.5cm)	3790
S	1	1	Marshall® Bladder Malleable Swivel Retractor, 3" x 6" (7.6cm x 15.2cm)	3792
T	2	2	Mayo Swivel Retractor, 2" x 2" (5cm x 5cm)	3883
U	2	2	Mayo Swivel Retractor, 2¾" x 2" (7cm x 5cm)	3335
V	--	--	Posterior Swivel Retractor, 1¾" x 3" (4.5cm x 7.6cm)	4975
W	--	--	Posterior Swivel Retractor, 1¾" x 4" (4.5cm x 10.2cm)	4976
X	1	1	Splanchnic Swivel Retractor, 1½" x 4" (3.8cm x 10.2cm)	3798



4849 White Bear Parkway, St. Paul, MN 55110 USA
 800-367-8657 • 651-287-4300 • fax: 651-287-4400
www.omni-tract.com

© 2006 Omni-Tract Surgical
 All Rights Reserved.

U.S. Patent Nos. 4,355,631; 4,617,916; 4,718,151; 4,930,932; 4,949,707;
 5,020,195; 5,242,240; 5,380,338; 5,400,772; 5,688,223; 5,727,899;
 5,792,046; 5,882,298; 5,899,627; 6,264,396; 6,315,718; 6,663,563,
 D343,235; Other Patents Pending. Omni-Flex and SideSnap are trademarks;
 Omni-Tract, Marshall and WISHBONE are registered trademarks of
 Minnesota Scientific, Inc.



5-06/40-6181/2500