



HlgH-Flo needles are designed as part of a complete subcutaneous infusion system, with options for every patient.



Please contact your healthcare provider if you have any questions regarding HlgH-Flo needles.

Name: \_\_\_\_\_

Telephone: \_\_\_\_\_

(write Pharmacy or Nurse information above)

**RMS resources (web links are case sensitive):**



Like us on **facebook**  
bit.ly/freedom60



**Top 10 FAQs**  
bit.ly/topFAQs



**Connections Blog Featuring Joanna Tierno**  
blog.rmsmedicalproducts.com



**Interactive Patient Information Guide**  
FREEDOM60<sup>®</sup>: joom.ag/VbXb  
FreedomEdge<sup>™</sup>: joom.ag/xffb



**6-minute FREEDOM60<sup>®</sup> SClg Training Video**  
bit.ly/adminvid1



**6-minute FreedomEdge<sup>™</sup> SClg Training Video**  
bit.ly/adminvid2



**Success with HlgH-Flo**  
rmsmedicalproducts.com/success

845-469-2042  
800-624-9600  
rmsmedicalproducts.com



HFL\_BR\_GetToKnow\_03-2046vG

The FREEDOM60<sup>®</sup> and FreedomEdge<sup>™</sup> Syringe Infusion Systems, RMS Precision Flow Rate Tubing<sup>™</sup> and HlgH-Flo Subcutaneous Safety Needle Sets<sup>™</sup> are registered trademarks of RMS Medical Products and are compliant with Medical Device Directive 93/42/EEC. RMS Medical Products is ISO 13485 certified. © 2016 RMS Medical Products; All Rights Reserved.



Get to *know*  
**HlgH-Flo**<sup>™</sup>  
Subcutaneous Safety Needle Sets

“It is easy to switch and I noticed improvements in my infusions right away.”

-Joanna T., SClg Patient




# Get to know the important components of HlgH-Flo

**1** Minimize tissue damage, scarring and local site reactions with the backcut needle tip design.

HlgH-Flo needle incision




other brand incision



backcut tip design

“The comfort of RMS HlgH-Flo needles is unparalleled... and I've tried them all!”

- Dave K., SCIG Patient & Specialty Pharmacy Sales Rep.



“The HlgH-Flo sub-q needles have a more consistent flow rate.”

- Gene B., RPh




custom luer design

**2**

Efficient & practical 20" tubing length:

- Less likely to tangle
- Faster flow rates



20"

“The HlgH-Flo needle tubing length is shorter than my previous product, and I like this... I feel more secure with the shorter tubing.”


- Maureen H., SCIG Patient

**Note:** 24" Extension set also available, for patients who need a longer length.

**3**

Needles stay in place after insertion, without tenting, due to the flexible, hinged butterfly wings.

“The butterfly wings on HlgH-Flo seem a bit broader than those of my previous product. This helps to secure the shorter 9mm needle in place so it doesn't dislodge during infusion.”



- Maureen H., SCIG Patient

**4**

Protect against accidental needle stick injuries with the safety closure feature.



“The higher-quality 3M Tegaderm™ dressing that RMS provides speaks to the way they care about their customers... a pleasant surprise.”

- Ted S., SCIG Patient

